

Features

- **Brushless Shade motor**
: Soft and long life cycle
- Various sample application
- Soft "Touch-On" mechanism
- Variable speed control
- Multi-function switch enables continuous or Touch-on control
- Safety functions
: Over-heat/load/current protection for safety of products and samples
- Made in Korea



Specifications

Model	VM1
Speed Range	0 ~ 3,300 rpm
Shaking Motion	Orbital (Circuit orbit 4.0mm)
Operation Mode	Touch or Continuous
Speed Control	Stepless
Motor	Brushless Shade motor
Head (Included)	Pop-off cup head and Φ 76mm Platform head
Max Load / Power	0.5kg/30W
Material	Aluminum die-casting with powder coating
Safety Device	Over-heat/load/current protection
Permissible Temp./RH %	5 ~ 40°C/80% RH
Dimensions / Weight	131x185x125mm/2.64kg
Packing Size / Gross Weight	190x250x215mm/2.94kg
Power Supply	AC120V, 60Hz or AC230V, 50/60Hz, 1P